

[KI] LLG-Magnetic stirrer RCT standard 2.0

The next generation of our top-selling IKA/LLG-Labware co-branded LLG-RCT Standard magnetic stirrer. New design, new features! Trusted and reliable quality!

- Digital controlled speed and temperature
- Clearly arranged LED display
- Digital display of set safety temperature limit
- Hot top indicator, hot surface warning to prevent burns
- Adjustable safety circuit of heating plate temperature
- Digital error code display
- Connection for external PT1000 temperature probe
- Aluminium hotplate for optimal heat transfer
- Enhanced safety due to improved heat control technology
- USB interface
- Powerful motor
- Improved magnetic coupling
- Supplied with housing protection cover
- Tested acc. to DIN EN IEC 61010-1

Scope of supply: Magnetic stirrer, temperature probe PT1000, screwdriver, protective cover, C16 mains cable (H11)



Max. stirring quantity (H ₂ O):	20 l
Speed range:	50 ... 1500 rpm
Heating plate:	Ø 135 mm
Temperature range:	RT + device self-heating ... 310 °C
Heat output:	600 W
Safety temperature range:	50 ... 360 °C
Dimensions (W x D x H):	160 x 270 x 85 mm
Weight:	2.4 kg
Power supply:	220 ... 230/115/100 V, 50/60 Hz
IP code:	IP 42
Warranty:	3 years

<i>Description</i>	<i>Stirring capacity</i>	<i>Min. speed</i>	<i>Max. speed</i>	<i>Plug type</i>	<i>PK</i>	<i>Cat. No.</i>	<i>Manufacturer Article Number</i>
	<i>l</i>	<i>rpm</i>	<i>rpm</i>				
LLG-Magnetic stirrer RCT standard 2.0	20	50	1500	EU	1	6.320400	6320400
LLG-Magnetic stirrer RCT standard 2.0	20	50	1500	UK	1	6.320402	6320402
LLG-Magnetic stirrer RCT standard 2.0	20	50	1500	CH	1	6.320401	6320401